

Custom Sling Design Script

Date: _____ Therapist: _____ Client: _____

1. CHOOSE A BASE SLING DESIGN (recommended) plus your preferred fixing.

Base Sling Design (loop fixing only)	
General Purpose	
General Purpose Head Support	
Care	
Care Head Support	
Junior	
Amputee	
Deluxe Leg	
Walking	

Base Sling Design	Loop	Clip	E-Clip	Dual (Clip)	Dual (E-Clip)
Pivot					
Pivot Care					
Standing Transfer					
Super Soft					
Other					

2. CHOOSE A BASE SLING SIZE. If you have selected "Other" please provide generic dimensions (A-D) plus any other relevant dimensions in the sketch box provided on the next page.

Base Sling Size							
Extra Small		Small		Medium		Large	
X Large		XX Large		Other			


3. CHOOSE A STANDARD MODIFICATION.

MOD 001 - add/change torso belt.

Torso Belt					
Buckle		Neoprene/Buckle Over		Custom Length	
Velcro		Neoprene/Velcro		Custom Width	
Comment					


MOD 002 - add/reduce head support

Head Support	
Add Head Support	cm
Reduce Head Support	cm




MOD 003 - add/reduce body length

Body Length	
Add Body Length	cm
Reduce Body Length	cm




MOD 004 - add/reduce body width

Body Width	
Add Body Width	cm
Reduce Body Width	cm




NB - Body width adjustment will not affect generic external profile of sling.


Leg Length	
Add Leg Length	cm
Reduce Leg Length	cm




Leg Width	
Add Leg Width	cm
Reduce Leg Width	cm




MOD 007 - Add Modesty Strap Modesty strap keeps legs of sling together without crossing.





MOD 008 - Care Sling Leg Strap Assists in support of upper thighs for clients with low tone using Care Slings.



MOD 009 - Stand-up Hemi Band Supports arm affected by hemiplegia in standing transfer.




MOD 010 - Stand-up Drop Strap Assists clients with limited weight bearing in standing transfer.



Sling Sketch and Dimensions

Custom modifications (please include in sketch)

A=



Disclaimer: While Custom Slings are covered by manufacturer's warranty, no responsibility is taken for their suitability and functionality.